KIT-GMP-G

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Cleanroom Gripping Solution

Gimatic electric grippers can operate in environments characterised by a high standard of cleanliness and hygiene by integrating the KIT-GMP-G. An adapter enables fastening the electric gripper to a robot interface (not included) and protecting its operation by means of a sterilisable silicone cover. The gripping fingers (not included) can be fastened by means of special hygienic stainless steel screws. The systems meets the most stringent GMP (Good Manufacturing Practice) cleanliness standards for Grade A and B.
Its main characteristics are:

- *Suitable for gripper MPPM1606.*
- *Compliant with GMP Grade A and B cleanliness standards.*
- *Cover made in Silicone Silpuran®Wacker FDA 21 CFR §177.2600.*
- *Compatible with hydrogen peroxide (VHP-H2O2) and UV cleaning procedures.*
- *Complete system with IP65 protection rating.*
- *The MPPM1606-KIT-GMP kit is certified to ISO14644-1 for ISO 2 cleanrooms, while the single gripper MPPM1606 is certified to ISO14644-1 for ISO 4 cleanrooms.*
- *Stainless steel hygienic screws (Novonox).*
- *Transparent surface that enables reading the LEDs of the internal sensors and the power supply status LEDs.*
- *Hygienic design surfaces preventing the build-up of bacteria.*
- *Easy fastening of gripping fingers with anti-rotation system.*

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







- *Fix robot interface to the robot by screws.*

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
- *Fix gripper to the gripper interface with screws and its centering sleeves forming the gripper subassembly.*
- *Maximum clamping torque of the screws = 1.5 Nm.*

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
- *Fix the gripper subassembly on the main interface.*

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
- *Fix internal jaws by screws and centering sleeves.*
- *At this point it is possible to add magnetic sensors on the groove to check jaws position. Cabling to the robot interface should be done in advance.*
- *Recommended clamping torque of the screws = 1.0 Nm.*
- *Maximum clamping torque of the screws = 1.5 Nm.*

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
- *Cover the system with the special silicone cap and place screw for fingers.*
- *Maximum clamping torque of the screws = 1.5 Nm.*

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
Users can create special gripping fingers for their specific application needs and fit them to the system by means of medical screws. The following figure shows an example of gripping fingers made with a two-screw fastening.
Fastening can be achieved with one or two screws depending on the shape.

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

Caution

- *Do not place the product in place where it may come in contact with foreign matter such as oil nor in an atmosphere of corrosive gas or flammable gases nor in place near inflammable material.*
- *Do not install the product near heating element.*
- *Do not perform wiring nor operate the product with wet hand.*
- *Do not transfer the product by holding only the cable.*
- *Do not attempt to force gripper jaws' movement by hand.*

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IPA Certification reference

The system has been declared suitable for use in hygienic areas by Fraunhofer IPA Institute as stated in report No. GI 1410-728